



The
Patent
Office

08 OCT 1999

PCT/GB99 / 03316



INVESTOR IN PEOPLE

The Patent Office
Concept House
Cardiff Road
Newport

South Wales
NR10 8QQ NOV 1999

WIPO PCT

I, the undersigned, being an officer duly authorised in accordance with Section 74(1) and (4) of the Deregulation & Contracting Out Act 1994, to sign and issue certificates on behalf of the Comptroller-General, hereby certify that annexed hereto is a true copy of the documents as originally filed in connection with the patent application identified therein.

In accordance with the Patents (Companies Re-registration) Rules 1982, if a company named in this certificate and any accompanying documents has re-registered under the Companies Act 1980 with the same name as that with which it was registered immediately before re-registration save for the substitution as, or inclusion as, the last part of the name of the words "public limited company" or their equivalents in Welsh, references to the name of the company in this certificate and any accompanying documents shall be treated as references to the name with which it is so re-registered.

In accordance with the rules, the words "public limited company" may be replaced by p.l.c., plc, P.L.C. or PLC.

Re-registration under the Companies Act does not constitute a new legal entity but merely subjects the company to certain additional company law rules.


Signed

Dated 10 November 1999

**PRIORITY
DOCUMENT**

SUBMITTED OR TRANSMITTED IN
COMPLIANCE WITH RULE 17.1(a) OR (b)

THIS PAGE BLANK (USPTO)



The Patent Office

1200198 E396465-3 002136
P01/7700 0.00 - 9822115.3

Patents Act 1977
(Rule 16)

Request for grant of a patent

(See the notes on the back of this form. You can also get an explanatory leaflet from the Patent Office to help you fill in this form)

The Patent Office

Cardiff Road
Newport
Gwent NP9 1RH

09 OCT 1998

1.	Your reference	SC/MH/N7444		
2.	Patent application number (The Patent Office will fill this part)	9822115.3		
3.	Full name, address and postcode of the or of each applicant (underline all surnames)	KING'S COLLEGE an Institution Incorporated by Royal Charter The Strand London WC2R 2LS		
	Patents ADP number (if you know it)	400185001		
	If the applicant is a corporate body, give the country/state of its incorporation	U.K.		
4.	Title of the invention	TREATMENT OF INFLAMMATORY DISEASE		
5.	Name of your agent (if you have one)	Williams, Powell & Associates		
	"Address for service" in the United Kingdom to which all correspondence should be sent (including the postcode)	4 St. Paul's Churchyard London EC4M 8AY		
	Patents ADP number (if you know it)	5830310001 ✓		
6.	If you are declaring priority from one or more earlier patent applications, give the country and the date of filing of the or of each of these earlier applications and (if you know it) the or each application number	Country	Priority application number (if you know it)	Date of filing (day / month / year)
7.	If this application is divided or otherwise derived from an earlier UK application, give the number and the filing date of the earlier application.	Number of earlier application		Date of filing (day / month / year)
8.	Is a statement of inventorship and of right to grant of a patent required in support of this request? (answer 'Yes if: a) any applicant named in part 3 is not an inventor, or b) there is an inventor who is not named as an applicant, or c) any named applicant is a corporate body. See note (d))	YES		

Patents Form 1/77

9. Enter the number of sheets for any of the following items you are filing with this form.
Do not count copies of the same document

Continuation sheets of this form

Description 23

Claim(s) 1

Abstract

Drawing(s)



10. If you are filing one of the following, state how many against each item.

Priority documents

Translations of priority documents

Statement of inventorship and right to grant of a patent (*Patents Form 7/77*)

Request for preliminary examination and search (*Patents Form 9/77*)

Request for substantive examination (*Patents Form 10/77*)

Any other documents
(please specify)

11. I/we request the grant of a patent on the basis of this application.

Signature

Date

William Anderson

08/10/98

12. Name and daytime telephone number of person to contact in the United Kingdom

Mr Lee Anderson 0171 329 4400

Warning

After an application for a patent has been filed, the Comptroller of the Patent Office will consider whether publication or communication of the invention should be prohibited or restricted under Section 22 of the Patents Act 1977. You will be informed if it is necessary to prohibit or restrict your invention in this way. Furthermore, if you live in the United Kingdom, Section 23 of the Patents Act 1977 stops you from applying for a patent abroad without first getting written permission from the Patent Office unless an application has been filed at least 6 weeks beforehand in the United Kingdom for a patent for the same invention and either no direction prohibiting publication or communication has been given, or any such direction has been revoked.

Notes

- If you need help to fill in this form or you have any questions, please contact the Patent Office on 0645 500505.*
- Write your answers in capital letters using black ink or you may type them.*
- If there is not enough space for all the relevant details on any part of this form, please continue on a separate sheet of paper and write "see continuation sheet" in the relevant part(s). Any continuation sheet should be attached to this form.*
- If you have answered 'Yes' Patents Form 7/77 will need to be filed.*
- Once you have filled in the form you must remember to sign and date it.*
- For details of the fee and ways to pay please contact the Patent Office.*

TREATMENT OF INFLAMMATORY DISEASE

This invention relates to inflammatory disease and more particularly to rheumatoid arthritis.

Rheumatoid arthritis (RA) is a chronic inflammatory disease of the synovial joints which leads to joint destruction, disability and early death. Although the cause of RA is presently unknown, it has been suggested that type II collagen, uniquely found in the articular cartilage, is a possible autoantigen for RA. It has recently been proposed that gp39, a 39KD glycoprotein, and peptides derived from it, are such autoantigens. However, the data supporting this hypothesis are limited and the role of gp39 therefore remains uncertain.

The present invention stems from a different approach based on a study of chondrocytes, the specialised cells of articular cartilage. We have isolated a protein from human chondrocytes and human chondrosarcoma cell lines which reacts with antibodies present in RA patients' sera and meets the accepted criteria for a putative autoantigen. This purified protein has been tested for proliferation of T cells and has been shown to selectively proliferate synovial T cells from patients with RA. This protein is the immunoglobulin heavy chain binding protein BiP(78KD). Correct recognition of an RA autoantigen leads to the development of prognostic and diagnostic tests for this disease and specific therapy. We have isolated and sequenced the DNA for this protein. We have also cloned and expressed this DNA. The amino acid and DNA sequences are novel and are shown in Appendices given hereafter.

The first part of the following description concerns the characterisation of such an autoantigen; the second the cloning, sequencing and expression of the protein; and the third part the demonstration of disease (rheumatoid arthritis) and tissue (synovial compartment) T cell specificity to the autoantigen.

Part 1: CHARACTERISATION OF AUTOANTIGEN

Chondrocytes/chondrosarcoma cells.

Chondrocytes were isolated from cartilage obtained during joint replacement. The cartilage was minced finely and digested with 1mg/ml collagenase. (Worthington). Following digestion the cells were centrifuged at 300g and resuspended in Dulbeccos minimal essential medium (DMEM)(Life Technologies, Paisley, UK) enriched with 10% foetal calf serum (FCS)(Harlan Sera-Lab, Loughborough, UK). Cell debris was washed off the adherent cells after 24hours and the cells allowed to expand until confluent. Cells were passaged using trypsin (0.25%) and split 1:3.

Chondrosarcoma cells (HTB94) (SW1353) were supplied by the American Type Culture Collection ATCC (Rockville, Maryland, USA) and by Dr J Block, Rush University, Chicago, USA (personal gift). These cells were cultured in DMEM with 10% FCS and split 1:3 after gentle trypsinisation (0.25% trypsin Life Technologies, Paisley, UK) when confluent.

Preparation of cell lysates.

Cells were scraped from the flask surface, homogenised and sonicated in the presence of proteinase inhibitors PMSF (2mM), leupeptin (200µg/ml) and aprotinin (50µg/ml) (Sigma, Poole, UK).

Sodium dodecyl sulphate (SDS) (Sigma, Poole, UK) was added to a final concentration of 1% and the proteins solubilised at room temperature for 1hour. Protein concentration was estimated by bicinchoninic acid assay using bovine serum albumin (BSA) as a standard protein (Sigma) and the cell lysate was used at 10µg/well equivalent.

Polyacrylamide gel electrophoresis (PAGE) and Western Blot.

The Mini Protean system (BioRad Laboratories, Hemel Hempstead, UK) was used to run the gels. 5, 7.5 or 10% SDS polyacrylamide denaturing gels 1.5mm thick were prepared (See Appendix 1). Gels were loaded with 10µg protein/well or the equivalent was loaded on preparative gels.

Electrophoresis was carried out at a constant 100V and broad range kaleidoscope markers (BioRad) were run in parallel with the cell lysates.

Following electrophoresis the proteins were blotted onto nitrocellulose at a constant 100V for 1 hour (See Appendix 1). Nitrocellulose was then blocked with 3%BSA (Sigma, Poole,UK) and left at 4⁰ C overnight. Preparative gel membranes were cut into 16 thin strips when necessary, each having an identical protein profile. The membranes were probed with patients sera (1/100 dilution) or specific monoclonal antibodies (at required concentrations) for 1 hour at room temperature and then washed x3 over 1 hour in TTBS (See Appendix 1). The secondary antibody, goat anti-human IgG (Fab²) horse radish peroxidase(HRP) conjugate (Sigma) was added at 1/1000 dilution and incubated for 1 hour at room temperature. The membranes were then washed x3 over 1 hour in TTBS. Enhanced chemiluminescence (Amersham) was used to develop the system and antigen-antibody.HRP complexes appeared as discrete bands on photographic film when developed.

Isolation of the putative autoantigen p78

The band of interest was seen in approximately 30% of the rheumatoid arthritis sera used to screen the cell lysates as previously described.

To isolate the protein the cell lysate was concentrated x 23, using a 30,000MW cut-off filter(Vivascience). This protein was then loaded on 5% and 7.5% gels in parallel . One lane on each was loaded at normal concentration while the two other lanes were overloaded with the protein. Kaleidoscope markers were loaded on either side of the test lanes. The gels were then run as previously described until the kaleidoscope markers showed that the 70,000 MW protein would be in the bottom third of the gel as close as possible to the run-off point. The gels were then blotted onto PVDF membrane (Immobilon P, Millipore) which was immediately placed in distilled water after transfer of the proteins was complete. The strip with normal loading was used for immunodetection of the protein band. The developed film of this immunodetection and Ponceau red staining of the overloaded strips was used to identify the band on the membrane which was then air-dried.

These strips were then taken to isolate and sequence the protein using matrix assisted laser desorption ionisation (MALDI) spectroscopy .

The electroblotted proteins were stained with Ponceau S (0.05% w/v aqueous methanol/0.1% acetic acid) using a rapid-staining protocol (1). The dried, stained proteins were then digested in situ with trypsin (Boehringer, modified) and the peptides extracted with 1:1 v/v formic acid:ethanol (2). One 0.2 µl aliquot (approximately 5% of the total digest) was sampled and directly analysed by matrix-assisted laser desorption ionisation (MALDI) time-of-flight mass spectrometry using a LaserMat 2000 mass spectrometer (Thermo Bioanalysis,UK) (3). A second 0.2µl aliquot was quantitatively esterified using 1% v/v thionyl chloride in methanol and also analysed by MALDI to provide acidic residue composition (4). Native and esterified peptide masses were then screened against the MOWSE peptide mass fingerprint database (5). The remaining digested peptides (90% of total digest) were then reacted with N-succinimidyl-2(3-pyridyl) acetate (SPA) in order to enhance b-ion abundance and facilitate sequence analysis by tandem mass spectrometry (6). Dried peptide fractions were treated with 7 µl 1% w/v N-succinimidyl-2(3-pyridyl) acetate in 0.5M HEPES (pH 7.8 with NaOH) containing 15% v/v acetonitrile for 20 min on ice. The reaction was terminated by 1 µl heptafluorobutyric acid (HFBA) and the solution immediately injected onto a capillary reverse-phase column (300 µm x 15cm) packed with POROS R2/H material (Perseptive Biosystems, MA) equilibrated with 2% v/v acetonitrile/0.05% v/v TFA running at 3 µl/min. The absorbed peptides were washed isocratically with 10% v/v acetonitrile/0.05% v/v TFA for 30 minutes at 3 µl/min to elute the excess reagent and HEPES buffer. The derivatised peptides were then eluted with a single step gradient to 75% v/v acetonitrile/0.1% v/v formic acid and collected in a single 4 µl fraction. The derivatised peptides were then sequenced by low-energy collision-activated dissociation (CAD) using a Finnigan MAT TSQ7000 fitted with a nanoelectrospray source (7,8). CAD was performed using 2.5 mTorr argon with collisional offset voltages between -18V and -28V. The product-ion spectra were collected with Q3 scanned at 500 amu/sec.

RESULTS:

Sequence data obtained from 7.5% gel (single band)

from GR78_human

Specific identifying peptides:

NQLTSNPENTVFDK 82-96
SDIDEIVLVGGSTR 353-366
TWNDPSVQQDIK 107-113

Identified human protein GR 78:
Kd glucose related protein precursor (GRP78)
Immunoglobulin heavy chain binding protein (BIP)

Part 2: CLONING, SEQUENCING AND EXPRESSION OF p78

1) mRNA isolation and PCR amplification of identified genes.

Human chondrocytes were isolated and cultured for three weeks as described. Poly(A) mRNA (1-2 μ g) was extracted with a Micro-Fastrack kit (Invitrogen) from a total of $1-2 \times 10^6$ cells. One microgram of the resulting mRNA was reverse transcribed into cDNA in a 20 μ l volume at 45°C for 1 hour using 1 μ l of Moloney murine leukemia virus reverse transcriptase (200u/ μ l); 5 X first strand buffer (Tris-HCl pH 8.3, 250mM; KCl 375 mM; MgCl 15mM); 0.1M DTT; oligo dT(12-18) 20ng/ μ l (Life Technologies); and dNTP mix 100mM (Amersham Pharmacia Biotech, Uppsala Sweden).

PCR was performed in a 50 μ l reaction volume under standard conditions (see below) using a Perkin Elmer Applied Biosystems thermal cycler PE2400. Primer sequences were derived from the GenBank database sequence corresponding to the Human gene for the immunoglobulin heavy chain binding protein, Bip (grp78), accession number X87949. Specific primers were synthesized to amplify the putative autoantigen gene from the chondrocyte cDNA. The resulting PCR product consisted of most of the grp78 coding region, bar the untranslated regions, signal sequences and the stop codon (nucleotide positions 279-2169 of the grp78 database sequence).

Primer sets for PCR were designed with integrated restriction sites to allow rapid subcloning of cDNA into the bacterial expression vector. The Forward primers encoded an NdeI site and the Reverse primers contained an XhoI restriction site:

Bip	Forward	primer	5'
TATACATATGGAGGAGGACAAGAAGGAGGACG 3' (32mer)			
Bip	Reverse	primer	5'
CCACCTCGAGTTCTGCTGTATCCTCTTCACCA 3' (32mer)			

After initial denaturation at 96°C for 2min the, PCR was performed for 28 cycles using a cycling profile of 94°C for 30s, 60°C for 30s and 72°C for 2min, with a final extension at 72°C for 7min. The PCR reaction generated a single specific Bip fragment of 1890bp.

2) Cloning of PCR generated fragments.

The restriction sites engineered into the forward and reverse primers used for the PCR reactions required flanking DNA for them to be recognised by their specific endonucleases (NdeI and XhoI). To provide this flanking DNA, the PCR generated fragment was cloned into a PCR cloning vector pCR2.1-TOPO (Invitrogen). The ligated plasmids were transformed into competent *E.coli* XL1-Blue (Stratagene) and plasmid DNA extracted using miniprep purification columns (Qiagen). The purified plasmid DNA for the clone was designated Bip-Topo. These DNA samples were stored at -20°C. The purified plasmid DNA for Bip-Topo was digested with NdeI and XhoI. The restricted fragments were separated by agarose gel electrophoresis and purified using the Qiagen DNA gel extraction kit.

3) Subcloning of restricted gene fragments into bacterial expression vector.

The purified fragment for the clone was ligated into the NdeI /XhoI pre-digested bacterial expression vector pET30a (Novagen). Ligation was performed at 12°C overnight in the presence of T4 ligase (20 units) and 1/10 vol of 10X ligase buffer (provided with the T4 ligase enzyme from Promega). The ligated plasmids were transformed into competent *E.coli* XL1-Blue (Stratagene) and screened by colony-PCR using Bip specific primers. Positive transformants carrying the required recombinant plasmids were purified and transformed into competent *E.coli* expression strain BL21-(DE3) (Invitrogen).

4) Sequencing of the 1890bp pET30::Bip subclone .

Sequencing of the 5' and 3' terminal ends of the pET30a::Bip clone confirmed that the recombinant DNA molecule was in-frame with the ATG start codon on the pET30 vector and that readthrough from this site continued through the Bip gene and ended with the 6X His residues and the stop codon located on the 3' arm of the expression vector.

Extensive DNA sequencing was performed using synthetic oligonucleotide primers spanning the entire length of the Bip subclone. Sequence analysis of the newly subcloned Bip gene fragment was performed by comparative alignment against the existing grp78 sequence from the database (accession number X87949). A number of differences between the two sequences were detected, both at the DNA and protein level (see Appendix 2). These areas of disagreement may either be a result of errors in the original DNA sequencing (of grp78) or they may indicate the presence of an additional related, but slightly different Bip gene in the genome.

All DNA sequencing was performed on an Applied Biosystems ABI 377 automated DNA sequencer using the dRhodamine dye terminator kit (Perkin Elmer- Applied Biosystems).

5) Expression of bacterial and purification recombinant proteins.

E.coli expression strain BL21-(DE3) containing the recombinant pET30a-Bip plasmid was grown at 37°C in LB medium containing kanamycin (50µg/ml). When the cells had reached an OD₆₀₀ of 0.6 units, isopropyl β-D-thiogalactopyranoside (IPTG) (1mM) was added to the medium to induce expression of the recombinant protein, driven by the IPTG-inducible promoter of the expression vector. To allow for maximal expression of the recombinant protein the culture was incubated for a further 4 hours at 37°C. The cells were pelleted by centrifugation and stored at -70°C.

For purification of the recombinant bacterial proteins the bacterial pellets were lysed in binding buffer (20mM NaPO₄, 500mM NaCl, 5mM Imidazole, 1mM PMSF, 1mg/ml Lysozyme, 5units/ml DNase, 0.1% Triton X-100, pH7.4). The lysate was cleared by centrifugation to remove insoluble matter and cell debris. The cleared lysate was passed

over a binding buffer-equilibrated chelating Hi-trap affinity column, with a bed volume of 5ml (Pharmacia). The non-specifically bound protein was washed from the column under stringent conditions using a series of three wash buffers. The primary wash was performed using 100ml of Binding buffer. This was followed by a high stringency low pH wash (20mM NaPO₄, 500mM NaCl, 0.1% Triton X-100, pH5.5) and an additional high stringency wash using 100ml of 20mM NaPO₄, 500mM NaCl, 0.1% Triton X-100, 50mM Imidazole, pH7.4.

The histidine-tagged recombinant proteins were eluted from the column by stripping with 50mM EDTA. Eluted proteins were dialysed overnight against 1x PBS to remove EDTA and Ni contaminants. The purified protein was concentrated and washed in sterile PBS using a 50000 Mw cutoff concentrator column (Millipore). The total amount of protein was determined by spectrophotometry using BSA as a standard with the bicinchoninic acid assay (according to manufacturers instructions).

The concentrated Bip recombinant protein was resuspended in PBS to a final concentration of 10mg/ml and stored at -70°C in 50 µl aliquots.

Part 3: T CELL SPECIFICITY IN RA TO p78

Lymphocyte proliferation and gamma-interferon production by peripheral blood and synovial fluid mononuclear cells to the putative autoantigen p78

Mononuclear cells were separated from heparinised peripheral blood or synovial fluid by density centrifugation using Lymphoprep (Nycomed-Amersham, Little Chalfont, UK). Cells were washed x3 with Hanks buffered saline solution then used at 2×10^5 /well. Cells were cultured in RPMI 1640 supplemented with 10% heat inactivated human serum. Proliferation assays were carried out using control stimulants, tuberculin protein purified derivative (PPD)(Central Veterinary Laboratories, Weybridge, UK) at optimal concentration (10µg/ml) and pokeweed mitogen (PWM)(Sigma, Poole, UK)(2µg/ml) and a range of concentrations for recombinant human p78 BiP (1-50µg/ml). Cells were cultured for six days at 37° C in 5% CO₂. 100µl of supernatants for IFN γ estimation were withdrawn on day 5 and fresh medium used to replenish wells. Tritiated thymidine (³HTdr) (0.2µCi/well) was added 24h prior to harvesting. Proliferation was determined by the uptake of

$^3\text{HTdr}$ estimated on a Dry Matrix beta counter (Canberra Packard) and expressed as counts per minute. Interferon γ was measured by ELISA. Paired monoclonal antibodies to Interferon γ (IFN γ) (Pharmingen, Oxford, UK) were used for capture ($1\mu\text{g/ml}$) and detection ($1\mu\text{g/ml}$) of IFN γ . A standard curve was run on each plate using recombinant human IFN γ over a range from $0.01\text{-}10\text{ngml}^{-1}$. Supernatants were used neat, 1/5, 1/10 and 1/100. The detection antibody was biotinylated, avidin-alkaline phosphatase conjugate was then linked to the biotin and developed with the substrate, 3,3',5,5'-tetramethylbenzidine (Sigma). The plates were read on a spectrophotometer at 450nm with a reference reading at 650nm.

1. T cell proliferation


T cell proliferative responses were determined for mononuclear cell preparations from paired peripheral blood and synovial fluid samples. Ten samples from patients with rheumatoid arthritis and ten from disease controls were examined. In addition six HLA-DR4 positive control samples were incorporated in the study. Responses were determined to an optimum concentration of p78 and control antigens. The control antigens included purified protein derivative of *Mycobacterium tuberculosis* (PPD), pokeweed mitogen (PWM) and β -galactosidase (produced by the same method as p78 as a control for bacterial contamination).

Results show a significant increase in the synovial fluid responses to p78 in the RA patients compared to peripheral blood (Stimulation index 4.8 ± 4.4 v 1.9 ± 1.5 ; $p=0.02$ Mann-Whitney test). These data suggest a specific role for this protein in the synovial compartment of patients with RA.

Production of gamma-interferon (γ -IFN) was also determined for these samples and paralleled lymphocyte proliferation.

The DR4 positive control data show no difference in peripheral blood response when compared to disease control samples. Similarly no significant response was seen to β -galactosidase in the RA synovial compartment, although responses were seen in control peripheral blood samples.

Proliferative response to p78



Proliferative responses to human p78 of mononuclear cells isolated from the peripheral blood and synovial fluid of patients with rheumatoid arthritis and other inflammatory joint diseases.

2. Use of tests for detection of antibodies to p78 in biological fluids or culture supernatants

Several techniques can be used, such as agglutination, Western blotting, and ELISA.

ELISA protocol for the detection of antibodies to p78 in sera

ELISA plates are coated half with p78 in bicarbonate buffer and half with buffer alone. The plates are blocked with bovine or human serum albumin (1% BSA or HAS) or other blocking protein before diluted sera are added to both the p78 coated half of the plate and also the non-coated half. After washing, biotin conjugated antibody to human immunoglobulin (polyclonal anti α , γ or μ) is added to the plate. Bound anti-immunoglobulin is detected with streptavidin-horse radish peroxidase and suitable substrate. Sera containing antibodies to p78 are determined spectrophotometrically. This test forms the basis of a diagnostic and/or prognostic test for rheumatoid arthritis.

3. Therapeutic application

There are several approaches to using p78 or derivatives for therapeutic purposes, including the following:

(a) Induction of mucosal tolerance.

Delivery of p78 autoantigen or peptides derived therefrom by mucosal routes, e.g. through the intestine or nasal mucosa, alters the immune response by downregulating disease activity leaving the patient's immune system otherwise intact. Alternatively p78 or p78 peptides can be delivered as DNA plasmids encoding them with an appropriate mammalian expression vector.

(b) Vaccination with TCR peptides

Peptides of the CDR3 region of the T cell receptor V α and V β chains can be synthesised and used as vaccines for delivery by intradermal or intramuscular injection. Plasmids encoding these peptides can be used in the same way.

(c) MHC blockade with native or altered peptides

The p78 peptides may be given parenterally or orally in appropriate cases either unmodified or modified by amino acid substitution and/or attachment of chemical groupings so as to block MHC and especially HLA-DR4 thereby leading to suppression of T cell activation and disease. P78 peptides either native or altered may be combined with soluble HLA-DR4 molecules and applied parenterally or orally.

(d) Induction of tolerance by plasmid DNA immunisation

Plasmids consisting of DNA coding for whole p78 protein or its peptides linked to a mammalian expression vector may be given by injection. DNA coding for human IL-10, IL-4, IL-11, or TGF-beta, incorporated singly or in any combination, may be used to deviate the immune response to p78 towards a TH2 mode so as to suppress disease.

In the therapeutic regimes indicated above the protein or derived peptide may be administered in appropriate compositions delivering amounts ranging from about 0.1 micrograms to about 1 gram or the equivalent in the case of plasmid or vaccine preparations.

Appendix 1

Methodology for gel electrophoresis

Acrylamide gel: 10%

6.075ml acrylamide (40%) (BDH, Poole, UK)
3.35ml methylenbisacrylamide (2%)(Pharmacia Biotech, Uppsala, Sweden)
6.25ml acrylamide gel buffer (see below)
9ml distilled water (produced within the laboratory)
250µl Ammonium persulphate (AMPERS)(0.025mg in 250ul of distilled water)(Sigma-Aldrich, Poole, UK)
25µl NNN'N'-Tetramethylethylenediamine (TEMED)(Sigma-Aldrich, Poole, UK)

Acrylamide gel buffer: pH 8.8

1.5M Tris (tris(hydroxymethyl)aminomethane)(Sigma-Aldrich, Poole, UK)
pH titration with concentrated hydrochloric acid
0.4% Sodium dodecyl sulphate (SDS)(BDH-Merck, Poole, UK)

Stacking gel:

1.2ml acrylamide (40%)(BDH-Merck, Poole, UK)
0.65ml bisacrylamide (2%)(Pharmacia Biotech, Uppsala, Sweden)
3.15ml stacking gel buffer (see below)
7.5ml distilled water (produced within the laboratory)
125µl ammonium persulphate (AMPERS) (0.025mg/250µl)(Sigma-Aldrich, Poole, UK)
12.5µl NNN'N'-Tetramethylethylenediamine (TEMED)(Sigma-Aldrich, Poole, UK)

Stacking gel buffer: pH 6.8:

0.5M Tris (tris(hydroxymethyl)aminomethane)(Sigma-Aldrich, Poole, UK)
pH titration with concentrated hydrochloric acid
0.4% Sodium dodecyl sulphate (SDS)(BDH-Merck, Poole, UK)

Loading buffer:

2ml glycerol
2ml 10% Sodium dodecyl sulphate (SDS)(BDH-Merck, Poole, UK)
0,25mg bromophenol blue
2.5ml stacking gel buffer 4-times concentrated (0.5M Tris; 0.4% SDS;
pH 6.8)
0.5ml 2-mercaptoethanol (Sigma-Aldrich, Poole, UK)

Electrophoresis/Running buffer:

3g/l Tris (tris(hydroxymethyl)aminomethane)(Sigma-Aldrich, Poole, UK)
14.4g/l glycine (BDH-Merck, Poole, UK)
1g/l Sodium dodecyl sulphate (SDS)(BDH-Merck, Poole, UK)

Transfer buffer:

3g/l Tris (tris(hydroxymethyl)aminomethane)(Sigma-Aldrich, Poole, UK)
14.4g/l glycine (BDH-Merck, Poole, UK)
1g/l Sodium dodecyl sulphate (SDS)(BDH-Merck, Poole, UK)
14.4g/l glycine (BDH-Merck, Poole, UK)
10% methanol (BDH-Merck, Poole, UK)

Tris Tween buffered saline (TTBS):

2.4g/l Tris (tris(hydroxymethyl)aminomethane)(Sigma-Aldrich, Poole, UK)
29g/l sodium chloride (BDH-Merck, Poole, UK)
500µl tween 20 (polyoxyethylene-sorbitan mono-laurate)(Sigma-Aldrich, Poole, UK)

3% Bovine serum albumin solution (BSA):

3g Albumin bovine fraction 5 (BDH-Merck, Poole, UK)
100ml TTBS (see above)

Appendix 2

Cloned pET30::Bip DNAa sequence aligned with grp78 DNA
(accession number X87949)

[illegible]Database RNA seq X87949
Cloned pET30BFP, seq

Database: JNA seq: X87949
Cloned: pET30BIP seq:

Database DNA seq. X87949

Database DNA seq X879499
Cloned pET30B1P seq

Database LNA seq X87949
Cloned pET30BIP seq

Database DNA seq X87949
Cloned pET30BIP seq

Database DNA seq X87949
Cloned pET30BIP seq

Database DNA seq X87949
Cloned pET30BIP seq

Database DNA seq X87945
Cloned pET30BIP seq

550	560	570	580	590	600
TCTGTGCAGCAGGACATCAAGTTCCTTGCCGCTTCAGGTTGGTTGA	AAAGAAAGAAACCTA	AAACCA			
TTCTGTGCAGCAGGACATCAAGTTCCTTGCCGCTTCAGGTTGGTTGA	AAAGAAAGAAACCTA	AAACCA			

Database DNA seq X87949
Cloned pET30BIP seq

610	620	630	640	650	660
TACATTCACAGTTTGATATTTGGAGGTGGGCAAAACAAGACATTTGCT	CTGAAGAAATTTCT				
TACATTCACAGTTTGATATTTGGAGGTGGGCAAAACAAGACATTTGCT	CTGAAGAAATTTCT				

Database DNA seq X87949
Cloned pET30BIP seq

670	680	690	700	710	720
GGCCATGGTTCCTCACTAAATAATGAAGAAGAAACCGCTGAGGCTTAT	TTGGGAAGAAAGTTAAC				
GGCCATGGTTCCTCACTAAATAATGAAGAAGAAACCGCTGAGGCTTAT	TTGGGAAGAAAGTTAAC				

Database DNA seq X87949
Cloned pET30BIP seq

730	740	750	760	770	780
CATGCAGTTGTACTGTACCAAGCCCTATTTTAATGATGCCCAACGCC	CAAGCAACCAAGAAC				
CATGCAGTTGTACTGTACCAAGCCCTATTTTAATGATGCCCAACGCC	CAAGCAACCAAGAAC				

Database DNA seq X87949
Cloned pET30BIP seq

790	800	810	820	830	840
GCTGGAAACCTATTCCTGGCCCTTAAATGTTATGAGGATCATCAACG	AGCCCTAACGGCAGCT				
GCTGGAAACCTATTCCTGGCCCTTAAATGTTATGAGGATCATCAACG	AGCCCTAACGGCAGCT				

Database DNA seq X87949
Cloned pET30BIP seq

850	860	870	880	890	900
ATTTGCTTATGGCCCTGGATTAAGAGGGAGGGGGAAGAAACATCCT	GGTGTGTGAACCTGGGT				
ATTTGCTTATGGCCCTGGATTAAGAGGGAGGGGGAAGAAACATCCT	GGTGTGTGAACCTGGGT				

Database DNA seq X87949
Cloned pET30BIP seq

910	920	930	940	950	960
GGCGGAACCTTCGATGTGTCTCTTCTCACCATTTGACCAATGGTG	CTTTCGAAGTTGTGGCC				
GGCGGAACCTTCGATGTGTCTCTTCTCACCATTTGACCAATGGTG	CTTTCGAAGTTGTGGCC				

Database DNA seq X87949
Cloned pET30BIP seq

970	980	990	1000	1010	1020
ACTAATGGAGATACCTCATCTGGGGTGGAGAAAGACTTTTGACCAAG	CGTGTCTCATGGAAACCTTC				
ACTAATGGAGATACCTCATCTGGGGTGGAGAAAGACTTTTGACCAAG	CGTGTCTCATGGAAACCTTC				

Database DNA seq X87949
Cloned pET30BIP seq

1030	1040	1050	1060	1070	1080
ATAAAGTGTACAAAGAGAGACGGGCAAAAGATGTCAAGGAAGGACAA	ATAGAGCTGTGCAG				
ATAAAGTGTGTACAAAGAGAGACGGGCAAAAGATGTCAAGGAAGGACAA	ATAGAGCTGTGCAG				

Database DNA seq X87949
Cloned pET30BIP seq

AACTCTTCAGATTCTTTTCTACACAGCTTCTGATAATCAACCAACTGTGTACCAATCAAGGCTCTAT Database DNA seq X87941
AAGTCTTCAGATTCTTTTCTACACAGCTTCTGATAATCAACCAACTGTGTACCAATCAAGGCTCTAT Cloned pET30BIP seq

FC

1.730

 \mathbb{J}, \mathbb{V}

1790

7.6

1850

1. V

1910

AG

1970

GGT

2030

66

2090

 $\bar{A} \bar{A}$

2150

2160	2170	2180	2190	2200	2210
A C A G C A G A A A A G A T G A G T T G T A G A C A C T G A T C T G C T A G T G C T G T A A T T G T A A A T					
Database DNA seq X87949					
Cloned pET30BIP seq					

2220	2230	2240	2250	2260	2270
A C T G G A C T C A G G A A C T T T T G T T A G G A A A A A T T G A A A G A A C T T A A G T C T C G A A T G T A A T T					
Database DNA seq X87949					
Cloned pET30BIP seq					

2300	2290	2300	2310	2320	2330
G G A A T C T T C A C C T C A G A G T G G A G T T G A A C T G C T A T A G C C T A A G C G G C T G T T A C T G C T T T					
Database DNA seq X87949					
Cloned pET30BIP seq					

2340	2350	2360	2370	2380	2390
T C A T T A G C A G T T G C T C A C A T G T C T T T G G G T G G G G G A G A A G A A A T T G G C C A T C T T A					
Database DNA seq X87949					
Cloned pET30BIP seq					

2400	2410	2420	2430	2440	2450
A A A A G C G G T A A A A A A C C T G G G T T A G G G T G T G T T C A C C T T C A A A A T G T T C T A T T T A A C					
Database DNA seq X87949					
Cloned pET30BIP seq					

2460	2470	2480	2490	2500	2510
A A C T G G G T C A T G T G C A T C T G G T G T A G G A A G T T T T C T A C C A T A A G T G A C A C C A A T A A T					
Database DNA seq X87949					
Cloned pET30BIP seq					

2520	2530	2540	2550
G T T T G T T A T T T A C A C T G G T C A A A A A A A A A A A A			
Database DNA seq X87949			
Cloned pET30BIP seq			

Decoration 'Decoratation #1': Shade (with solid black) residues that differ from Database DNA seq X87949.

Appendix 3

Expressed pET30::Bip amino acid sequence aligned with predicted
protein from grp78 DNA sequence (accession X87949)

	10	20	30	40	50	60	
1	L	S	V	A	A	M	L
2	L	L	L	S	A	R	A
3	E	E	E	D	K	K	E
4	D	K	K	E	D	V	G
5	T	V	G	I	D	L	G
6	T	T	Y	S	C	V	G
7	F	K	N	G	R	V	E
8	I	I	A	N	D	Q	G
9	N	R	P	E	T	I	30BIP
10	expressed	protein					
11	E	E	D	K	K	E	D
12	V	G	I	D	L	G	T
13	T	Y	S	C	V	G	F
14	K	N	G	R	V	E	I
15	I	A	N	D	Q	G	N
16	R	P	E	T	I	30BIP	expressed
17	protein						
18							
19							
20							
21							
22							
23							
24							
25							
26							
27							
28							
29							
30							
31							
32							
33							
34							
35							
36							
37							
38							
39							
40							
41							
42							
43							
44							
45							
46							
47							
48							
49							
50							
51							
52							
53							
54							
55							
56							
57							
58							
59							
60							
61							
62							
63							
64							
65							
66							
67							
68							
69							
70							
71							
72							
73							
74							
75							
76							
77							
78							

Alignment Report of Unfilled, using Clustal method with PAM250 residue weight table.
 Tuesday, September 8, 1998 12:58 pm

	560	560	570	580	590	
140 HVND A E K F A E E D K K L K E R I D T R N E L E S Y A Y S L K N Q I G D K E K L G G K L S S E D K E T M E K A V E E						grp78 X87949 protein
523 HVND A E K F A E E D K K L K E R I D T R N E L E S Y A Y S L K N Q I G D K E K L G G K L S S E D K E T M E K A V E E						PEP30BIP expressed protein
600	610	620	630			
600 K I E W L E S H Q D A D I E D F K A K K K E L E E I V Q P I I S K L Y G S A G P P P T G E E D T A E K D E L						grp78 X87949 protein
583 K I E W L E S H Q D A D I E D F K A K K K E L E E I V Q P I I S K L Y G						PEP30BIP expressed protein

Insertion 'D' : Shade (with solid black) residues that differ from grp78 X87949 protein.

CLAIMS

1. The use of immunoglobulin heavy chain binding protein BiP (78KD), or a peptide derived therefrom, as a marker for the presence of rheumatoid arthritis.
2. The use of protein BiP(78KD) according to claim 1, in which the protein has the amino-acid sequence as shown in Appendix 3.
3. The protein BiP(78KD) having the amino-acid sequence shown in Appendix 3.
4. A prognostic or diagnostic test for RA in which the test reagent for testing body fluid is immunoglobulin heavy chain binding protein BiP (78KD), or a peptide derived therefrom.
5. A test according to claim 4, being an ELISA assay.
6. Immunoglobulin heavy chain binding protein BiP (78KD), or a peptide derived therefrom, for use in therapy.
7. Immunoglobulin heavy chain binding protein BiP (78KD), or a peptide derived therefrom, for use in rheumatoid arthritis therapy.
8. A DNA sequence coding for Immunoglobulin heavy chain binding protein BiP (78KD), or for a peptide derived therefrom, for use in therapy.
9. A DNA sequence coding for Immunoglobulin heavy chain binding protein BiP (78KD), or for a peptide derived therefrom, for use in rheumatoid arthritis therapy.
10. A DNA sequence according to claim 9, containing the nucleotide sequence having the lower lines of sequence information shown in Appendix 2.
11. A recombinant vector containing a DNA sequence according to claim 9 or 10.